



CERTIFICATE OF ANALYSIS

**Product: DMEM-LOW GLICOSE – 1 X Concentrated
With Sodium Bicarbonate (3.7 g/L)**

Lot # 180.317BT

Catalog # BR30002.05

Manufacture Date: Mar/13/2017

Expiration Date: Mar/13/2018

Test	Specification	Units	Results
Appearance	Clear, transparent		Clear, transparent
pH w/ NaHCO ₃	7.2 +/- 0.2		7.21
Bacteria and yeast	Negative		Negative
Endotoxins	≤ 10.0	EU/mL	≤ 1,0
<i>Mycoplasma</i>	Negative	qPCR	Negative
Osmolality w/NaHCO ₃	305 +/- 5%	mOsm/Kg H ₂ O	311

Responsável Técnico
LGC Biotecnologia Ltda.

Rua Pasadena, 235 – Pq. Empresarial San José
Cotia – SP / Brazil - 06715-864
Phone-Fax: 55 11 46148070

Rua Pasadena, 235 – Pq. Empresarial San José
Cotia – SP / Brazil - 06715-864
Phone-Fax: 55 11 46148070