
CERTIFICATE OF ANALYSIS

**Product: DMEM-LOW GLICOSE – 1 X Concentrated
With Sodium Bicarbonate (3.7 g/L)**


Lot # 141.117B

Catalog # BR30002.05

Manufacture Date: Nov/09/2017

Expiration Date: Nov/09/2018

Test	Specification	Units	Results
Appearance	Clear, transparent		Clear, transparent
pH w/ NaHCO ₃	7.2 +/- 0.2		7.32
Bacteria and yeast	Negative		Negative
Endotoxins	≤ 10.0	EU/mL	≤ 2,0
<i>Mycoplasma</i>	Negative		Negative
Osmolality w/NaHCO ₃	305 +/- 5%	mOsm/Kg H ₂ O	310



Responsável Técnico
LGC Biotecnologia Ltda.