
CERTIFICATE OF ANALYSIS


Product: MEM-ALPHA MODIFICADO
With ribo and desoxiribonucleosides

Lot # 040.717TJB

Catalog # BR30007.05

Manufacture Date: Jun/29/2017
Expiration Date: Jun/29/2018

Test	Specification	Units	Results
Appearance	Clear, transparent		Clear, transparent
pH w/ NaHCO ₃	7.2 +/- 0.2		7.23
Bacteria and yeast	Negative		Negative
Endotoxins	≤ 10.0	EU/mL	≤ 2,0
<i>Mycoplasma</i>	Negative		Negative
Osmolality w/NaHCO ₃	255 +/- 5%	mOsm/Kg H ₂ O	265



Responsável Técnico
LGC Biotecnologia Ltda.